Magnetic Frame

Product Name: Magnetic Frame

Package Size: 1 set

Application: For separation of magnetic beads and solution.

Description: These products are suitable for the magnetic separation of magnetic beads in 1.5 ml centrifuge tube, 15 ml centrifuge tube, 200 µl PCR plate, 8-tube strip or 200 µl PCR tubes. The average magnetic separation time is 10 sec. It can be used for the rapid separation of nucleic acids, proteins, cells and other substances.

Product Components:

These products consist of the following components:

| Reagent Name | OSE-MF-01 | OSE-MF-02 | OSE-MF-03 | OSE-MF-04 |
|-----------------------------------|-----------|-----------|-----------|-----------|
| Magnetic Frame | 1 set | - | - | - |
| Demagnetization Stand | 1 set | - | - | - |
| 1.5 ml Centrifuge Tube Holder | 1 set | - | - | - |
| 15 ml Centrifuge Tube Holder | 1 set | 1 set | - | - |
| 15 ml Centrifuge Tube Base | 1 set | - | - | - |
| 200 µl PCR Plate Magnetic Frame | - | - | 1 set | - |
| 8-Tube Strip (0.2 ml) Magic Frame | - | - | - | 1 set |
| Handbook | 1 | 1 | 1 | 1 |

Precautions: These products contain very strong permanent magnets, please be caution about the following notes.

1. People wearing a pacemaker or any other medical magnetized implant should not use this product unless advised by a health professional.
2. Try to avoid contact with metal substances and other magnets.
3. Keep the product away from credit cards, all magnetic products and recording equipment.
4. Keep the product away from precision instruments, watches, electronic equipment, displays and monitors.
5. These products have strong magnetism and can attract ferromagnetic materials such as iron, nickel, cobalt and other metals, so it should be handled carefully.
6. The applicable operating temperature of the product is -40°C ~ 80°C.

Cleaning and Disinfection

These products are not suitable for high temperature and high pressure sterilization. If cleaning is needed, spray and/or wipe with 70% ethanol, 1% sodium hypochlorite solution or 0.1 M HCl solution. Do not submerge in aqueous solution for a long time.

Manufacturer

Manufacturer name:
TIANGEN BIOTECH (BEIJING) CO., LTD.

Registered and production address:
Building 5, No. 86, Shuangying West Road, Science and Technology Park, Changping District, Beijing

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