

## Taq DNA Polymerase

Cat. no. 4992761/4992762/4992764/4992763/4992765

Storage: -30~-15°C for two years.

Concentration: 2.5 and 5 U/μl

Product size

Product Components	Taq DNA Polymerase	10×Taq Buffer	10×Taq Buffer (Mg <sup>2+</sup> free)	MgCl <sub>2</sub> (25 mM)
4992761	250 U (2.5 U/μl)	1.8 ml	/	/
4992762	500 U (2.5 U/μl)	1.8 ml	/	/
4992764	500 U (5 U/μl)	1.8 ml	/	/
4992763	500 U (2.5 U/μl)	/	1.8 ml	1.8 ml
4992765	500 U (5 U/μl)	/	1.8 ml	1.8 ml

TIANGEN BIOTECH (BEIJING) CO., LTD.  
[HTTP://WWW.TIANGEN.COM/EN](http://www.tiangen.com/en)

The product is used for research only, neither intended for the diagnosis, or treatment of a disease, nor for the food, or cosmetics etc.

### Introduction

Taq DNA Polymerase is a recombinant 94 kDa DNA polymerase expressed in an *E.coli* strain that carries the cloned *Thermu aquaticus* DNA Polymerase gene. It possesses both 5'-3' polymerase and exonuclease activity, and has no detectable 3'-5' exonuclease activity. The extension rate of Taq DNA Polymerase is 1-2 kb/min in PCR. In addition, it has a 3' adenylation activity. Thus, the PCR products can be used directly in TA-cloning procedures.

### Unit Definition

One unit of Taq DNA Polymerase is defined as the amount that incorporates 10 nmol of dNTPs into acid-insoluble material within 30 min at 74°C with activated salmon sperm DNA as the template/primer.

### Quality Control

Purity up to standard by SDS-PAGE assay. No foreign nuclease activity is detected. Single copy gene in human genome could be amplified effectively. No significant activity change after storing at room temperature (15-30°C) for one week.

### Storage Buffer

20 mM Tris-HCl (pH 8.0), 0.1 mM EDTA, 1 mM DTT, 100 mM KCl, Stabilizers, 50% Glycerol.

### 10×Taq Buffer

200 mM Tris-HCl (pH9.0), 200 mM KCl, 100 mM (NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub>, 15 mM MgCl<sub>2</sub>, other components.

- Two types of 10× Taq buffer can be chosen: Mg<sup>2+</sup> plus and Mg<sup>2+</sup> free.
- Mg<sup>2+</sup> free buffer is supplied with separate 25 mM MgCl<sub>2</sub> Solution.
- Unless specifically requested, Mg<sup>2+</sup> plus buffer will be supplied as regular component.

### Applications

PCR amplification of DNA fragments, DNA labeling, primer extension, DNA sequencing and addition of extra adenines at the 3'-end of PCR products. PCR products with A-tailing at 3' are suitable for TA cloning.

### Example

**Note: The following example is only for reference, user must set up optimal reaction system according to different reaction conditions such as different templates or primers etc.**

1. For 50 μl PCR reaction system: 1 kb fragment of human genomic DNA was amplified (If use different reaction system, please proportionally increase or decrease the amount of reaction components referring to this system).

Template	< 1 μg
Primer 1(10 μM)	1 μl
Primer 2(10 μM)	1 μl
10×Taq Buffer	5 μl
dNTP Mixture(2.5 mM)	4 μl
Taq (2.5 U/μl)	0.5-1 μl
ddH <sub>2</sub> O	up to 50 μl

2. PCR cycle set-up:

94°C 3 min	} 30 cycles
94°C 30 sec	
55°C 30 sec	
72°C 1 min	
72°C 5 min	

3. Result detection: Load 5 μl PCR products to agrose gel for detection.