

TGuide Smart Tissue/Cell DNA Kit

(Prefilled 96-Deepwell plate)

For genomic DNA purification from tissue and cells.

TECHNICAL MANUAL

Cat. no. GDP602-E

Note: To use the TGuide Smart Tissue/Cell DNA Kit, you must have the TGuide Smart Tissue/Cell DNA (program no. DP602) installed on the TGuide S16/S32 Nucleic Acid Extractor.



www.tiangen.com/en

This product is for scientific research use only.
Do not use in medicine, clinical treatment, food
or cosmetic

Table Contents

<i>Kit Contents</i>	<i>1</i>
<i>Tissue/Cell DNA reagent composition</i>	<i>1</i>
<i>Storage condition</i>	<i>1</i>
<i>Product</i>	<i>2</i>
<i>Features</i>	<i>2</i>
<i>Notes</i>	<i>2</i>
<i>Operational steps</i>	<i>3</i>
<i>1.Prefilled 96-Deepwell plate</i>	<i>4</i>
<i>2.Sample processing</i>	<i>4</i>
<i>3.Operation steps of TGuide S16 Nucleic Acid Extractor</i>	<i>5</i>
<i>Appendix</i>	<i>6</i>
<i>1.Program</i>	<i>6</i>
<i>2.Related Products</i>	<i>6</i>

TGuide Smart Tissue/Cell DNA Kit

Cat. no. GDP602-E

Kit Contents

Contents	GDP602-E (96 preps)
Buffer GHA	50 ml
Tissue/Cell DNA Reagents	6 plates
Proteinase K	2×1 ml
RNase A (100 mg/ml)	400 µl
Buffer TB	15 ml
TGuide Smart Tip Comb	12 pcs

Tissue/Cell DNA reagents composition

Column 1/7	Column 2/8	Column 3/9	Column 4/10	Column 5/11	Column 6/12
Buffer GHLP	Buffer GDAP	Buffer GDAP	Buffer PWDP	None	MagAttract Suspension GSP1
600 µl	900 µl	900 µl	900 µl		520 µl

Storage condition

The kit can be stored under dry conditions at room temperature (15~30°C) for 12 months. If the solution precipitates, it can be preheated in a water bath at 37°C for 10 min to dissolve the precipitation, without affecting the effect.

Product

This kit adopts magnetic beads and a unique buffer system to separate and purify high-quality genomic DNA from various animal tissues and cells. The uniquely embedded magnetic beads have a strong affinity for nucleic acid under certain conditions. When the conditions are changed, the magnetic beads can release the absorbed nucleic acid to rapidly separate and purify the nucleic acid.

It can be used to perfectly fit with TGuide S16 Nucleic Acid Extractor. Through absorption, transfer and release of magnetic beads by the special magnetic bar, magnetic beads and nucleic acid can be transferred to improve the degree of automation. The whole process is safe and convenient, and the extracted genomic DNA fragments are large, with high purity and reliable quality.

The DNA purified by this kit is suitable for a range of common downstream applications including digestion, PCR, library construction, Southern hybridization, and other experiments.


Features

- **Simple and fast:** Ultra-pure genomic DNA can be obtained by running TGuide S16 for 50 minutes.
- **Wide use:** It is applicable to all kinds of animal tissue and cells.
- **Ultra-pure:** The obtained DNA has high purity and can be directly used in PCR, digestion, hybridization and other molecular biological experiments.

Note

1. Repeated freezing and thawing samples should be avoided, otherwise the extracted DNA fragments will be small and the total yield will decrease.
2. If there is precipitation in the Buffer GHA, it can be dissolved in a 37°C water bath and used after shaking well.

Operational steps



1. Prefilled 96-Deepwell plate

1.1 Take out a prefilled 96-Deepwell plate and invert it to re-suspend the magnetic beads; Gently shake to concentrate the reagent and magnetic beads to the bottom of the plate. Before use, remove sealing film carefully to avoid liquid spatter or spills.

1.2 Add proper volume (60~100 μ l) of buffer TB to Column 5/11 of the plate.

2. Sample processing

2.1 Tissue

Take 10~50 mg of animal tissue, use a liquid nitrogen or high-throughput tissue grinding homogenizer (TGrinder H24R tissue homogenizer, self-prepared, TIANGEN: OSE-TH-02) to adequately grind the tissue, or cut tissue into small pieces as much as possible, then add 400 μ l Buffer GHA and 20 μ l Proteinase K. Proceed to 3.1.

1) Samples with tissue blocks visible to naked eyes are recommended to be digested at 65°C for 30 min until completely digested;

2) Samples with sufficient homogenization do not need the above digestion process.

3) The rat tail samples should be digested overnight at 56°C.

Notes: After sample digestion, if there are tissue fragments, it is recommended to centrifuge at 12,000 rpm for 1 min to remove residual impurities. To remove RNA, add 4 μ l RNaseA (100 mg/mL), shake it for 15 sec, and incubate at room temperature for 5 min.

2.2 Cell

1) Processing methods of different cell samples:

(1) Suspension cell: Determine the number of cells collected (the collected number should not be more than 1×10^7) and centrifugate at $300 \times g$ for 5 min. Then collect cells into a centrifuge tube and carefully remove all supernatant of the culture medium. Wash the cells with PBS solution, and then suck out the PBS solution as much as possible. Add 100 μ l PBS to the cells, and completely re-suspend the cells.

(2) Adherent cell: Determine the number of cells and remove the culture medium. Wash cells with the PBS solution, and suck out the PBS solution. Then add the PBS solution containing 0.10~0.25% trypsin to cells for digestion. When the cells are released from the wall of the vessel, add a medium containing serum to inactivate the trypsin. Transfer the cell solution to an RNase-free centrifuge tube and centrifuge it at $300 \times g$ for 5 min. Collect cell pellets and carefully remove all supernatant. Add 100 μ l PBS to the cells, and completely re-suspend the cells.

Note: When collecting cells, it is important to remove all cell culture medium; otherwise, it will lead to incomplete digestion.

- 2) Add 200 μ l Buffer GHA and 20 μ l Proteinase K to the collected cell pellets, and completely re-suspend the cells.


Note. To remove RNA, add 4 μ l RNaseA (100 mg/mL), shake it for 15 sec, and incubate at room temperature for 5 min.


- 3) Proceed to 3.1.

3. Operation steps of TGuide S16 Nucleic Acid Extractor

- 3.1 Add 300 μ l solution after processing the above sample in Column 1/7 of the plate. Put cartridges on the reagent tank bracket of TGuide S16 Nucleic Acid Extractor.

- 3.2 Place the reagent plate on the base in the TGuide S16 Nucleic Acid Extractor. Insert the Tip Combs into slots to ensure that they are well connected and firmed.



- 3.3 If you use the TGuide S16 Nucleic Acid Extractor, select the corresponding program DP602 file on the touch screen, click the icon  in the lower right corner of the screen, or click the "RUN" button at the bottom of the screen to start the experiment.

- 3.4 At the end of the automated extraction process, take the DNA out of Column 5/11 of the plate and store it under appropriate conditions. 96-Deepwell plate and tip comb are for single use only.

Detection of DNA concentration and purity

The size of the obtained genomic DNA fragment is affected by the sample storage time and shear force during operation. The concentration and purity can be detected by agarose gel electrophoresis and UV spectrophotometer. Ideally, the DNA should absorb at most at OD₂₆₀, where an OD₂₆₀ value of 1 corresponds to approximately 50 μ g/ml double strand DNA and 40 μ g/ml single strand DNA. The OD₂₆₀/OD₂₈₀ ratio should be 1.7~1.9.

Appendix

1. Program

The extraction process of S16 provided for DP602 is shown in the following table:

Step	Hole site	Step name	Mix time (min)	Mix speed	Dry time (min)	Volume (μl)	Temp. (°C)	Segments	Every time (sec)	Magnetization time(sec)	Cycle	Magnet speed (mm/s)
1	1	Lysis	2	8	0	900	--	1	0	0	0	--
2	6	Collect beads	0.5	8	0	500	--	5	3	0	2	2.5
3	1	Bind	10	8	0	900	--	5	4	0	2	2.5
4	2	Wash 1	5	7	0	900	--	5	3	0	2	2.5
5	3	Wash 2	5	7	0	900	--	5	3	0	2	2.5
6	4	Wash 3	5	7	5	900	--	5	3	0	2	2.5
7	5	Elution	10	7	0	100	75	5	5	0	2	2.5
8	6	Discard	0.5	5	0	500	--	1	0	0	0	--

2. Related Products

Instrument and Accessories

Product name	Packing Size	Cat.No
TGuide S16 Nucleic Acid Extractor	1 set	OSE-S16
TGuide Smart Magnetic Tip Comb	200 pieces/box	OSE-TGA-S03
TGuide Single Sample Tank Bracket	5 pieces/box	OSE-TGA-S32

TGuide Smart Reagent Kits

Product name	Cartridge	Cat.No
TGuide Smart Magnetic Plant DNA Kit	48 preps	GDP607-DE
TGuide Smart Soil/Stool DNA Kit	48 preps	GDP612-DE
TGuide Smart Magnetic Plant RNA Kit	48 preps	GDP662-DE
TGuide Smart DNA Purification Kit	48 preps	GDP642-DE
TGuide Smart Blood/Cell/Tissue RNA Kit	48 preps	GDP661-DE
TGuide Smart Blood Genomic DNA Kit	48 preps	GDP601-DE
TGuide Smart Viral DNA/RNA Kit	48 preps	GDP604-DE
TGuide Smart Universal DNA Kit	48 preps	GDP605-DE
TGuide Smart Magnetic Tissue/Cell DNA Kit	96 preps	GDP602-E
TGuide Smart Blood/Cell Tissue RNA Kit	96 preps	GDP661-E
TGuide Smart Universal DNA Kit	96 preps	GDP605-E