



(DP304) TIANamp Genomic DNA Kit

—Animal Tissue

TIANGEN BIOTECH(BEIJING)CO., LTD

Ver. No. 20170328

WWW.TIANGEN.COM


Experiment Preparation

1. Rat hepatic 10-30 mg; Grinder
2. Pipette and matched sterile tips (10 μ l, 200 μ l, 1 ml); 1.5 ml centrifuge tubes
3. PBS solution, 96-100% ethanol
4. Vortex oscillator; Dry bath/water bath; Centrifuge



Experiment Preparation - Kit Preparation


Please add 96-100% ethanol in Buffer PW and GD before use according to the volume indicated on the label of the bottle.



Step 1




Add PBS (about 50 μ l) to the sample, and break into cell suspension with a tissue grinder




Centrifuge at 10,000 rpm ($\sim 11,500\times g$) for 1 min, and remove supernatant. Add 200 μ l Buffer GA and mix thoroughly by vortex.

Step 2



Add 200 μ l Buffer GB and 20 μ l Proteinase K, and mix well upside down.

Step 3



Place it at 70°C for 10 min, and the solution should be clear. Briefly centrifuge to remove water drops of the cap and inner wall.

Step 4



Add 200 μ l 96-100% ethanol, mix thoroughly by vortex for 15 sec, there may be flocculent precipitate at this time. Briefly centrifuge to remove water drops of the cap and inner wall.

Step 5



Transfer the solution and flocculent precipitate from the previous step to Spin Column CB3.



Centrifuge at 12,000 rpm ($\sim 13,400\times g$) for 30 sec, discard the waste liquid in the collection tube and place the Spin Column CB3 into the collection tube.

Step 6



Centrifuge at 12,000 rpm ($\sim 13,400\times g$) for 30 sec, discard the waste liquid in the collection tube and place the Spin Column CB3 into the collection tube.

Add 500 μ l Buffer GD to Spin Column CB3

Step 7



Add 600 μ l Buffer PW to Spin Column CB3 (ensure 96-100% ethanol has been added before use).



Centrifuge at 12,000 rpm (\sim 13,400 \times g) for 30 sec, discard waste liquid in the collection tube and place the Spin Column CB3 into the collection tube.

Step 8 Repeat step 7.

Step 9




Centrifuge at 12,000 rpm ($\sim 13,400 \times g$) for 2 min, and discard the waste liquid.



Place the Spin Column CB3 at room temperature for 2 minutes to completely dry Buffer PW in the membrane.

Note: Ethanol residues in Buffer PW can inhibit subsequent enzymatic reactions (restriction enzyme digestion, PCR, etc.) experiments.

Step 10



Transfer the Spin Column CB3 into a 1.5 ml centrifuge tube, and add 50-200 μ l Buffer TE to the middle of the adsorption membrane. Place at room temperature for 2 min and centrifuge at 12,000 rpm ($\sim 13,400\times g$) for 2 min to collect the solution into the centrifuge tube.