



# (DP316) TIANamp Micro DNA Kit

## ——Serum/Plasma


**TIANGEN BIOTECH(BEIJING)CO., LTD**

Ver. No. 20170330

[WWW.TIANGEN.COM](http://WWW.TIANGEN.COM)

# Experiment Preparation

1. Serum/plasma samples (this experiment takes serum as an example)
2. 96-100% ethanol
3. Pipette and matched sterile tips (10  $\mu$ l, 200  $\mu$ l, 1 ml); 1.5 ml centrifuge tubes
4. Vortex oscillator; Dry bath/water bath; Centrifuge



# Experiment Preparation-Kit Preparation


Please add 96-100% ethanol in Buffer PW and GD before use according to the volume indicated on the label of the bottle.



## Preparation of Carrier RNA stock solution

When using carrier RNA for the first time, please dissolve Carrier RNA (310  $\mu\text{g}$ ) in 310  $\mu\text{l}$  RNase-Free ddH<sub>2</sub>O, then aliquot the solution and store at -20°C. The concentration of the solution is 1  $\mu\text{g}/\mu\text{l}$ . Avoid repeated freezing and thawing for more than 3 times.

# Step 1



Take 100  $\mu$ l – 200  $\mu$ l serum/plasma into 2 ml centrifuge tube. If less than 100  $\mu$ l, add Buffer GA to 100  $\mu$ l final volume.

## Step 2



Add 20  $\mu$ l Proteinase K, and mix well by vortex. Add 200  $\mu$ l Buffer GB, and gently mix well upside down.

Incubate at 56°C for 10 min, and shake the sample during the period. Briefly centrifuge to remove water drops of the cap and inner wall.


(If serum/plasma volume <50  $\mu$ l, add 1  $\mu$ l of 1  $\mu$ g/ $\mu$ l Carrier RNA storage solution)

## Step 3



Add 200  $\mu$ l ethanol (96-100%). If the room temperature is higher than 25°C, please precool ethanol on ice. Gently mix the sample well upside down, and place at room temperature for 5 minutes. Briefly centrifuge to remove water drops of the cap and inner wall.

## Step 4



Transfer the solution from the previous step to Spin Column CR2. (place CR2 into the collection tube).

Centrifuge at 12,000 rpm (~13,400×g) for 30 sec, discard the waste liquid in the collection tube and place the Spin Column CR2 into the collection tube.

## Step 5



Add 500  $\mu$ l Buffer GD to Spin Column CR2 (make sure 96-100% ethanol has been added before use)

Centrifuge at 12,000 rpm ( $\sim 13,400\times g$ ) for 30 sec, discard the waste liquid in the collection tube and place the Spin Column CR2 into the collection tube.



## Step 6



Add 600  $\mu$ l Buffer PW to Spin Column CR2 (make sure 96-100% ethanol has been added before use).


Centrifuge at 12,000 rpm ( $\sim 13,400\times g$ ) for 30 sec, discard the waste liquid in the collection tube and place the Spin Column CR2 into the collection tube.

## Step 7 Repeat Step 6.

## Step 8




Centrifuge at 12,000 rpm ( $\sim 13,400\times g$ ) for 2 min, and discard the waste liquid.



Place the Spin Column CR2 at room temperature for 2 minutes to completely dry Buffer PW in the membrane.

## Step 9



Transfer the Spin Column CR2 into a 1.5 ml centrifuge tube, and add 20-50  $\mu$ l Buffer TB to the middle of the adsorption membrane. Place at room temperature for 2 min and centrifuge at 12,000 rpm ( $\sim 13,400 \times g$ ) for 2 min to collect the solution into the centrifuge tube.