



(DP321) DNAquick Plant System -Plant

TIANGEN BIOTECH(BEIJING)CO., LTD

Ver. No. 20170328

WWW.TIANGEN.COM

Experiment Preparation


1. Plant Leaves
2. Mortar; Liquid nitrogen
3. Pipette and matched sterile tips (10 μ l, 200 μ l, 1 ml); 1.5 ml centrifuge tubes
4. Isopropanol; 70% ethanol; Clean absorbent paper
5. Vortex oscillator; Dry bath/water bath; Centrifuge



Step 1



Add liquid nitrogen to 100 mg fresh plant tissue or 20 mg dry weight tissue, and fully grind the tissues.




Add 400 μ l Buffer FP1 and 6 μ l RNase A (10 mg/ml). After vortex for 1 min, place at room temperature for 10 min.

Step 2




Add 130 μ l Buffer FP2, mix well by vortex for 1min.

Step 3



Centrifuge at 12,000 rpm (~13,400 g) for 5 min, and transfer the supernatant to a new centrifuge tube.

Step 4



Add $0.7 \times$ volume of isopropanol to the supernatant and mix well. There may be flocculent genome DNA at this time. (for example, add 350 μ l isopropanol to 500 μ l supernatant),



Centrifuge at 12,000 rpm (\sim 13,400 g) for 2 min, discard supernatant and retain precipitation.

Step 5



Centrifuge at 12,000 rpm (~13,400 g) for 2 min, and discard the supernatant.

Add 600 μ l buffer 70% ethanol, and mix well by vortex for 5 sec.

Step 6 重复操作步骤5。

Step 6



Open cover and place on clean absorbent paper at room temperature for 5-10 minutes to completely dry the ethanol.

Note: Ethanol residues can inhibit subsequent enzymatic reactions (restriction enzyme digestion, PCR, etc.) experiments.

Step 7



Add appropriate amount of Buffer TE, and incubate at 65°C for 10-60 min to dissolve DNA. Mix well upside down for several times during the period.