



(DP330) TIANquick FFPE DNA Kit

——FFPE


TIANGEN BIOTECH(BEIJING)CO., LTD

Ver. No. 20170412

WWW.TIANGEN.COM


Experiment Preparation

1. Paraffin section or block
2. 96-100% ethanol
3. Pipette and matched sterile tips (10 μ l, 200 μ l, 1 ml); 1.5 ml centrifuge tubes
4. Vortex oscillator; Dry bath/water bath; Centrifuge



Experiment Preparation-Kit Preparation

Please add 96-100% ethanol in Buffer PW and GD before use according to the volume indicated on the label of the bottle.



Step 1



Take 5-8 pieces of paraffin sections (5-10 μm , $1 \times 1 \text{ cm}^2$).

Note: If the surface of the sample is exposed to air, the first 2~3 scraped pieces shall be discarded.

Step 2



Place the sample in a 1.5 ml sterile centrifuge tube, add 500 μ l Buffer GL, and then add 50 μ l Buffer GP. Vortex violently for 10 sec.

Step 3




Incubate at 98 °C for 30 min, mix well upside down for 3 times during the period, until the sample is completely dissolved.

Step 4



Centrifuge at 12,000 rpm (~13,400 g) for 5 min

Step 5



Pipette the water phase clear liquid of the middle layer along the tube wall with 200 μ l tip into the new centrifugal tube

(The upper layer is a mixture of paraffin and protein, while the lower layer is a little impurity precipitation. If the layering is not complete, the centrifugation time can be prolonged until the upper layer is separated from the aqueous phase clear liquid.).

Step 6



Add twice volume of 96-100% ethanol

Mix well and let it stand for 3 minutes

Step 7



Transfer the solution and flocculent precipitate from the previous step to Spin Column CR2. Centrifuge at 8,000 rpm ($\sim 6,000\times g$) at room temperature for 2 min, discard the waste liquid in the collection tube and place the Spin Column CR2 into the collection tube.

Step 8



Add 500 μ l Buffer GD to Spin Column CR2.

Centrifuge at 8,000 rpm ($\sim 6,000\times g$) at room temperature for 1 min, discard the waste liquid and place the Spin Column CR2 into the collection tube.

Step 9




Add 600 μ l Buffer PW to Spin Column CR2.



Centrifuge at 8,000 rpm ($\sim 6,000\times g$) at room temperature for 1 min, discard the waste liquid and place the Spin Column CR2 into the collection tube.

Step 10 Repeat step 9.

Step 11




Centrifuge at 12,000 rpm ($\sim 13,400\times g$) for 2 min, and discard the waste liquid.

Place the Spin Column CR2 at room temperature for 2-5 minutes to completely dry Buffer PW in the membrane.

Note: Ethanol residues in Buffer PW can inhibit subsequent enzymatic reactions (restriction enzyme digestion, PCR, etc.) experiments.

Step 12



Transfer the Spin Column CR2 into a 1.5 ml centrifuge tube, and add 65°C preheated 30-100 μ l Buffer TE or ddH₂O to the middle of the adsorption membrane. Placed at room temperature for 2-5 min and centrifuge at 12,000 rpm (\sim 13,400 \times g) for 2 min to collect the solution into the centrifuge tube. The purified DNA can be stored at -20°C.