



(DP334) TIANamp Blood Spots DNA Kit

——Blood Spots


TIANGEN BIOTECH(BEIJING)CO., LTD

Ver. No. 20170328

WWW.TIANGEN.COM


Experiment Preparation

1. Dry blood spot; Puncher (3 mm)
2. 96-100% ethanol
3. Pipette and matched sterile tips (10 μ l, 200 μ l, 1 ml); 1.5 ml centrifuge tubes
4. Vortex oscillator; Dry bath/water bath; Centrifuge




Experiment Preparation - Kit Preparation

Please add 96-100% ethanol in Buffer PW and GD before use according to the volume indicated on the label of the bottle.




Step 1



Add 3 pieces of dried blood spot samples with a diameter of 3 mm into a 1.5 ml centrifugal tube.

Step2




Add 200 μ l of Buffer GA.

Step 3



Add 20 μ l
Proteinase K




Vortex for 10
sec to mix



Incubate at 56°C for 1 hr, during which vortex
for 10 sec every 10 min.



Step 4



Add 200 μ l of Buffer GB, and mix thoroughly by vortex for 10 sec




Incubate at 70°C for 10 min, during which vortex for 10 sec every 3 min. Briefly centrifuge to remove water drops of the cap and inner wall after the incubation.

Step 5



Add 100 μ l 96-100% ethanol. If room temperature is higher than 25°C, please prechill ethanol on ice. Gently mix the sample upside down, and incubate at room temperature for 5 min. Briefly centrifuge to remove water drops of the cap and inner wall.

Step 6



Transfer the solution from the previous step to Spin Column CR2 (place CR2 in the collection tube).

Centrifuge at 12,000 rpm ($\sim 13,400\times g$) for 30 sec, discard the waste liquid in the collection tube and place the Spin Column CB2 into the collection tube.

Step 7



Centrifuge at 12,000 rpm (~13,400×g) for 30 sec, discard the waste liquid in the collection tube and place the Spin Column CR2 into the collection tube.

Add 500 µl Buffer GD to Spin Column CR2 (ensure that 96-100% ethanol has been added before use).

Step 8



Centrifuge at 12,000 rpm (~13,400×g) for 30 sec, discard waste liquid in the collection tube and place the Spin Column CR2 into the collection tube.

Add 700 µl Buffer PW to Spin Column CR2 (ensure 96-100% ethanol has been added before use).

Step 9



Centrifuge at 12,000 rpm ($\sim 13,400 \times g$) for 30 sec, discard waste liquid in the collection tube.

Add 500 μ l Buffer PW to Spin Column CR2 (ensure 96-100% ethanol has been added before use).

Step 10




Place the Spin Column CR2 in the collection tube. Centrifuge at 12,000 rpm ($\sim 13,400 \times g$) for 2 min, and discard the waste liquid.



Place the Spin Column CR2 at room temperature for 2-5 minutes to completely dry Buffer PW in the membrane.

Note: Ethanol residues in Buffer PW can inhibit subsequent enzymatic reactions (restriction enzyme digestion, PCR, etc.) experiments.

Step 10



Transfer the Spin Column CR2 into a 1.5 ml centrifuge tube, and add 20-50 μ l Buffer TB to the middle of the adsorption membrane. Place at room temperature for 2 min and centrifuge at 12,000 rpm ($\sim 13,400\times g$) for 2 min to collect the solution into the centrifuge tube.