




# TGuide Smart Processed Food DNA Kit

(Prefilled single sample cartridge)

*For extraction of genomic DNA from processed foods, including soybean powder, potato chips, starch, honey, and dried meat.*

## TECHNICAL MANUAL

Cat. no. GDP680-DE

**Note:** To use the TGuide Smart Processed Food DNA Kit, you must have the TGuide Smart Processed Food DNA Kit (program no. DP680) installed on the TGuide S16/S32/S32 Pro/S96 Dex Nucleic Acid Extractor.



[www.tiangen.com/en](http://www.tiangen.com/en)

For research use only. Not for use in diagnostic procedures.

## Table Contents

|                                                                                                          |                  |
|----------------------------------------------------------------------------------------------------------|------------------|
| <b><i>Kit Contents</i></b>                                                                               | <b><i>1</i></b>  |
| <b><i>Processed Food DNA Reagents Composition</i></b>                                                    | <b><i>1</i></b>  |
| <b><i>Storage</i></b>                                                                                    | <b><i>1</i></b>  |
| <b><i>Introduction</i></b>                                                                               | <b><i>2</i></b>  |
| <b><i>Features</i></b>                                                                                   | <b><i>2</i></b>  |
| <b><i>Notes</i></b>                                                                                      | <b><i>2</i></b>  |
| <b><i>Protocol</i></b>                                                                                   | <b><i>3</i></b>  |
| <b><i>I. Preparation of Processed Food DNA Extraction Reagents</i></b>                                   | <b><i>3</i></b>  |
| <b><i>II. Sample Pretreatment</i></b>                                                                    | <b><i>3</i></b>  |
| <b><i>III. Operation Steps for TGuide S16/S32/S32 Pro/S96 Dex<br/>        Nucleic Acid Extractor</i></b> | <b><i>6</i></b>  |
| <b><i>Appendix</i></b>                                                                                   | <b><i>10</i></b> |

# TGuide Smart Processed Food DNA Kit

Cat. no. GDP680-DE

## Kit Contents

|             | Contents                    | GDP680-DE<br>(48 preps) |
|-------------|-----------------------------|-------------------------|
| GDP680-DH   | Buffer GHA                  | 40 ml                   |
|             | Buffer GHL                  | 20 ml                   |
|             | NKY tablets                 | 2 tablets               |
|             | Processed Food DNA Reagents | 48 pcs                  |
|             | Proteinase K                | 1 ml                    |
|             | RNaseA (100 mg/ml)          | 200 µl                  |
|             | Buffer TBC                  | 15 ml                   |
| OSE-TGA-S36 | TGuide Smart Tip Comb       | 12 pcs                  |

**Note: GDP680-DH and OSE-TGA-S36 are shipped and packaged separately.**

## Processed Food DNA Reagents Composition

| Well 1    | Well 2    | Well 3    | Well 4     | Well 5 | Well 6                          |
|-----------|-----------|-----------|------------|--------|---------------------------------|
| Buffer IC | Buffer W3 | Buffer W3 | Buffer BP3 | None   | MagAttract<br>Suspension BEUSP5 |
| 450 µl    | 900 µl    | 900 µl    | 900 µl     |        | 915 µl                          |

## Storage

Store at room temperature (15-30°C) in a dry environment. If precipitates form in any buffer, incubate the solution in a 37°C water bath for 10 minutes to fully dissolve precipitates.

## Introduction

Processed foods have complex compositions, containing not only multiple raw material components but also various food additives (such as salt, sugar, oils, and pigments). Furthermore, processing techniques such as frying, deep-frying, boiling, and roasting can cause varying degrees of DNA degradation in the raw materials. Conventional nucleic acid extraction methods are unable to simultaneously achieve high yield, purity, and scalable extraction.

This kit employs a magnetic bead and buffer system with unique separation properties to isolate high-quality genomic DNA from a wide range of processed foods, including soybean powder, potato chips, starch, honey, and dried meat. The specially embedded magnetic beads exhibit strong affinity for nucleic acids under specific conditions and release them upon condition change, enabling rapid purification. The procedure is safe and convenient, yielding DNA fragments with high yield, excellent purity, and reliable quality. The kit is compatible with TIANGEN's TGuide S16/S32/S32 Pro/S96 Dex systems for automated extraction.

DNA purified with this kit is suitable for various downstream assays, including PCR, real-time PCR, etc.

## Features

**Simple and fast:** simple pretreatment; high-quality genomic DNA can be obtained within 1 hour after starting the instrument run.

**Safe and non-toxic:** no phenol/chloroform or other hazardous reagents required.

**Broad applicability:** suitable for genomic DNA extraction from various processed food samples, including soybean powder, potato chips, starch, honey, dried meat, etc.

**High quality:** the extracted DNA is of high quality and can be directly used in PCR, real-time PCR, and other assays.

## Notes Please read these notes carefully before using this kit.

1. If precipitates appear in Buffer GHA, redissolve them in a 37°C water bath and mix well before use.
2. The NKY tablets provided in the kit are an optional addition. For samples with high pigment or inhibitor content, such as brown sugar, strawberry jam, tomato ketchup, hawthorn jelly, instant noodle sauce packets, and lotus root starch, it is recommended to add the NKY tablets.

Two methods are available for using the NKY tablets; choose either one:

Method 1: First add the NKY tablet to Buffer GHA, with 1 tablet per 20 ml of solution. Vortex for 2 min to disperse — it will appear as a white, insoluble turbid suspension. Mix thoroughly before adding the lysis buffer. Slight unevenness when dispensing into centrifuge tubes does not affect the results.

Method 2: Crush the NKY tablet and add approximately a mung bean-sized particle to each sample, then add Buffer GHA and proceed with the operating steps.

3. Pretreatment methods vary for different sample types. Please read the manual carefully before the experiment and prepare the required reagents.

## Protocol

### I. Preparation of Processed Food DNA Extraction Reagents

- 1.1 Remove the single-cartridge prefilled reagents from the kit. Invert several times to resuspend the magnetic beads. Flick the strip gently to collect all reagents and beads at the bottom of the wells. Carefully remove the sealing film before use. Avoid vibration to prevent liquid splashing.
- 1.2 Add appropriate volume (60-100  $\mu$ l) of elution Buffer TBC to well 5 before starting the run.

### II. Sample Pretreatment

Refer to the table below for the recommended sample amount of different sample types:

| Sample Classification                                               | Recommended Sample Amount | Sample Name                                                                                                                                                                |
|---------------------------------------------------------------------|---------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| soy milk powder, milk powder, flour, and starch                     | 50-100 mg                 | soybean meal, soy milk powder, milk powder, sesame paste, starch, lotus root starch, vermicelli, rice noodles                                                              |
|                                                                     |                           | biscuits, bread, fried dough twist, senbei, snow rice cake, instant noodle cake, crispy rice crust                                                                         |
| meat, seafood, processed fruits and vegetables, sugar, and pet food | 50-100 mg                 | diced pork, dried pork slice, dried small shrimp, deep-sea dried fish, fish tofu, shredded squid, duck blood curd, sausage                                                 |
|                                                                     |                           | raisins, dried figs, instant noodle vegetable packet, pickled vegetables, potato chips, corn flakes, dried banana chips, oatmeal, konjac snack, dried tofu, cheese, butter |
|                                                                     |                           | marshmallow, crispy candy, White Rabbit creamy candy, chocolate, milk tablet                                                                                               |
|                                                                     |                           | cat food, cat treats                                                                                                                                                       |
| meat sauce, tomato ketchup, and yogurt                              | 100-200 mg                | salad dressing, mushroom sauce, instant noodle sauce packet, fermented bean curd, bean paste                                                                               |
| soy milk, milk, fruit juice, and honey                              | 1-10 ml                   | soy milk, milk, yogurt, concentrated fruit juice, honey                                                                                                                    |

Refer to the following procedure for sample pretreatment:

### 1. Soy Milk Powder, Milk Powder, Flour, and Starch

- 1) Place 50-100 mg of sample into a grinding tube (user-provided; TIANGEN, Cat. no. OSE-TH-B05) or a steel bead grinding tube (user-provided; TIANGEN, Cat. no. OSE-TH-B04). Add 700  $\mu$ l of Buffer GHA, 300  $\mu$ l of Buffer GHL, and 20  $\mu$ l of Proteinase K. Use a thermostatic shaking metal bath (user-provided; TIANGEN, Cat. no. OSE-DB-03) and digest at 65°C while shaking for 30 min.

**Note:**

- a. **The amount of Buffer GHA can be adjusted according to the sample's water absorbency.**
  - b. **For samples that are difficult to homogenize, such as vermicelli and brown rice, grind and mix using a TGrinder H24 Tissue Homogenizer (user-provided; TIANGEN, Cat. no. OSE-TH-01) before adding Buffer GHA (vibrate at 6 m/s for 30 sec, with a 30 sec interval, for 3-9 cycles).**
  - c. **For starch and vermicelli samples, the addition of lysis buffer may result in a jelly-like consistency. It is recommended to add 4  $\mu$ l of TR Amylase (user-provided; TIANGEN, Cat. no. GRT404-01) to aid digestion.**
- 2) After cooling to room temperature, add 4  $\mu$ l RNaseA (100 mg/ml), vortex to mix, and incubate at room temperature for 10 min.
  - 3) Centrifuge at 13,400 rpm ( $\sim$ 17,000 $\times$ g) for 5 min and add no more than 500  $\mu$ l of the digested solution to Well 1 of the single sample reagent cartridge for automated extraction.

### 2. Meat, Seafood, Processed Fruits and Vegetables, Sugar, and Pet Food

- 1) Flash-freeze 50-100 mg of sample in liquid nitrogen and place it into a grinding tube (user-provided; TIANGEN, Cat. no. OSE-TH-C01 or OSE-TH-C02) or a steel bead grinding tube (user-provided; TIANGEN, Cat. no. OSE-TH-B04). Use a TGrinder H24 Tissue Homogenizer (user-provided; TIANGEN, Cat. no. OSE-TH-01) to grind and mix (vibrate at 6 m/s for 30 sec, with a 30 sec interval, for 3-9 cycles). Then add 700  $\mu$ l of Buffer GHA, 300  $\mu$ l of Buffer GHL, and 20  $\mu$ l of Proteinase K, and place in a thermostatic shaking metal bath (user-provided; TIANGEN, Cat. no. OSE-DB-03) and digest at 65°C while shaking for 30 min.

**Note: For samples that are difficult to homogenize, such as certain dried meat products, increase the homogenization cycles until the sample is ground to a powder, or extend the 65°C digestion time until no visible tissue clumps remain.**

- 2) After cooling to room temperature, add 4  $\mu$ l RNaseA (100 mg/ml), vortex to mix, and incubate at room temperature for 10 min.
- 3) Centrifuge at 13,400 rpm ( $\sim$ 17,000 $\times$ g) for 5 min and add no more than 500  $\mu$ l of the digested solution to Well 1 of the single sample reagent cartridge for

automated extraction.

### 3. Meat Sauce, Tomato Ketchup, and Yogurt

- 1) Take 100-200 mg of sample, add 700  $\mu$ l of Buffer GHA, 300  $\mu$ l of Buffer GHL, and 20  $\mu$ l of Proteinase K, then place in a thermostatic shaking metal bath (user-provided; TIANGEN, Cat. no. OSE-DB-03) and digest at 65°C while shaking for 30 min. If visible clumps are present, transfer to a grinding tube (user-provided; TIANGEN, Cat. no. OSE-TH-C02 or OSE-TH-C03) or a steel bead grinding tube (user-provided; TIANGEN, Cat. no. OSE-TH-B04), use a TGrinder H24 Tissue Homogenizer (user-provided; TIANGEN, Cat. no. OSE-TH-01) to grind and mix (vibrate at 6 m/s for 30 sec, with a 30 sec interval, for 6-9 cycles), then add 700  $\mu$ l of Buffer GHA, 300  $\mu$ l of Buffer GHL, and 20  $\mu$ l of Proteinase K, and place in a thermostatic shaking metal bath and digest at 65°C while shaking for 30 min. If no visible clumps are present, the sample may also be added directly to a grinding tube (user-provided; TIANGEN, Cat. no. OSE-TH-B05), which facilitates nucleic acid extraction from microorganisms present in the sample.

**Note: The water content of different samples may vary slightly. Centrifugation may be performed as needed. For example, a 200 mg yogurt sample can be centrifuged at 5,000 $\times$ g for 10 min, and the supernatant discarded before proceeding with subsequent processing.**

- 2) After cooling to room temperature, add 4  $\mu$ l RNaseA (100 mg/ml), vortex to mix, and incubate at room temperature for 10 min.
- 3) Centrifuge at 13,400 rpm (~17,000 $\times$ g) for 5 min and add no more than 500  $\mu$ l of the digested solution to Well 1 of the single sample reagent cartridge for automated extraction.

### 4. Soy Milk, Milk, Fruit Juice, and Honey


- 1) Take 1-10 ml of sample, centrifuge at 12,000 rpm for 2 min, retain the pellet and 200  $\mu$ l of supernatant. Add 300  $\mu$ l of Buffer GHA, 300  $\mu$ l of Buffer GHL, and 20  $\mu$ l of Proteinase K, then place in a thermostatic shaking metal bath (user-provided; TIANGEN, Cat. no. OSE-DB-03) and digest at 65°C while shaking for 30 min.

**Note: Due to differences in concentration and content among samples, the volume may be adjusted accordingly. For example, concentrated fruit juice may require approximately 1 ml, while relatively clear juice may require about 10 ml.**

- 2) After cooling to room temperature, add 4  $\mu$ l RNaseA (100 mg/ml), vortex to mix, and incubate at room temperature for 10 min.
- 3) Centrifuge at 13,400 rpm (~17,000 $\times$ g) for 5 min and add no more than 500  $\mu$ l of the digested solution to Well 1 of the single sample reagent cartridge for automated extraction.

### III. Operation Steps for TGuide S16/S32/S32 Pro/S96 Dex Nucleic Acid Extractor

1. Add no more than 500 µl of the digested solution to Well 1 of the single sample reagent cartridge and place the single sample reagent cartridge onto the dedicated single sample reagent cartridge holder.
2. Mount the single sample reagent cartridge holder onto the base platform inside TGuide S16/S32/S32 Pro/S96 Dex instrument; insert tip combs into corresponding comb rack slots and fully lock the buckles in place.



3. Select the corresponding program and press the Start button to initiate the extraction protocol.

**The automated extraction program for processed food DNA is shown in the table below:**

Table 1 TGuide S16 Automated Nucleic Acid Extractor – Automated Extraction Program for Processed Food DNA

| Step | Slot | Name           | Mixing Time (min) | Mixing speed | Air-Drying Time (min) | Volume (μl) | Temp. (°C) | Magnetic Attraction Segments | Magnetic Attraction Time Per Segment (sec) | Liquid Surface Magnetic Attraction Time (sec) | Cycle Number | Magnetic Attraction Speed (mm/s) |
|------|------|----------------|-------------------|--------------|-----------------------|-------------|------------|------------------------------|--------------------------------------------|-----------------------------------------------|--------------|----------------------------------|
| 1    | 6    | Transfer Beads | 0                 | 8            | --                    | 930         | --         | 5                            | 5                                          | 0                                             | 2            | 2.5                              |
| 2    | 2    | Store Beads    | 1                 | 8            | --                    | 900         | --         | 1                            | 0                                          | 0                                             | 1            | 2.5                              |
| 3    | 1    | Binding        | 1                 | 8            | --                    | 900         | --         | 1                            | 0                                          | 0                                             | 1            | 2.5                              |
| 4    | 2    | Transfer Beads | 0                 | 8            | --                    | 930         | --         | 5                            | 5                                          | 0                                             | 2            | 2.5                              |
| 5    | 1    | Binding        | 3                 | 8            | --                    | 900         | --         | 1                            | 0                                          | 0                                             | 1            | 2.5                              |
| 6    | 1    | Binding        | 5                 | 5            | --                    | 900         | --         | 3                            | 20                                         | 10                                            | 1            | 2.5                              |
| 7    | 2    | Wash 1         | 0.5               | 8            | --                    | 100         | --         | 1                            | 0                                          | 0                                             | 1            | 2.5                              |
| 8    | 2    | Wash 2         | 4                 | 8            | --                    | 900         | --         | 5                            | 5                                          | 0                                             | 2            | 2.5                              |
| 9    | 3    | Wash 3         | 0.5               | 8            | --                    | 100         | --         | 1                            | 0                                          | 0                                             | 1            | 2.5                              |
| 10   | 3    | Wash 4         | 4                 | 8            | --                    | 900         | --         | 5                            | 5                                          | 0                                             | 2            | 2.5                              |
| 11   | 4    | Wash 5         | 3                 | 8            | --                    | 900         | --         | 5                            | 5                                          | 0                                             | 2            | 2.5                              |
| 12   | 6    | Wash 6         | 3                 | 6            | 8                     | 900         | --         | 3                            | 20                                         | 0                                             | 1            | 2.5                              |
| 13   | 5    | Elution        | 8                 | 8            | --                    | 120         | 75         | 5                            | 12                                         | 0                                             | 3            | 2.5                              |
| 14   | 6    | Discard Beads  | 0.5               | 8            | --                    | 300         | --         | 1                            | 0                                          | 0                                             | 1            | 2.5                              |

**Table 2: TGuide S32/S32 Pro Automated Nucleic Acid Extractor — Automated Extraction Program for Processed Food DNA**

| Step | Slot | Name           | Waiting time (min) | Mixing time (min) | Magnetic Attraction time (min) | Mixing Speed | Volume (μl) | Temperature (°C) | Cycle Mode |
|------|------|----------------|--------------------|-------------------|--------------------------------|--------------|-------------|------------------|------------|
| 1    | 6    | Transfer Beads | 0                  | 0                 | 1                              | Fast         | 930         | --               | Cycle      |
| 2    | 2    | Store Beads    | 0                  | 1                 | 0                              | Fast         | 900         | --               | --         |
| 3    | 1    | Binding        | 0                  | 1                 | 0                              | Fast         | 900         | --               | --         |
| 4    | 2    | Transfer Beads | 0                  | 0                 | 1                              | Fast         | 930         | --               | Cycle      |
| 5    | 1    | Binding        | 0                  | 3                 | 1                              | Fast         | 900         | --               | Cycle      |
| 6    | 1    | Binding        | 0                  | 5                 | 2                              | Slow         | 900         | --               | Cycle      |
| 7    | 2    | Wash 1         | 0                  | 0.5               | 1                              | Fast         | 100         | --               | Cycle      |
| 8    | 2    | Wash 2         | 0                  | 4                 | 1                              | Fast         | 900         | --               | Cycle      |
| 9    | 3    | Wash 3         | 0                  | 0.5               | 1                              | Fast         | 100         | --               | Cycle      |
| 10   | 3    | Wash 4         | 0                  | 4                 | 1                              | Fast         | 900         | --               | Cycle      |
| 11   | 4    | Wash 5         | 0                  | 3                 | 1                              | Fast         | 900         | --               | Cycle      |
| 12   | 6    | Wash 6         | 0                  | 3                 | 1                              | Slow         | 900         | --               | Cycle      |
| 13   | 5    | Elution        | 8                  | 8                 | 3                              | Fast         | 120         | 75               | Cycle      |
| 14   | 6    | Discard Beads  | 0                  | 0.5               | 0                              | Fast         | 300         | --               | --         |

**Table 3: TGuide S96 Dex Automated Nucleic Acid Extractor – Automated Extraction Program for Processed Food DNA**

| Step | Well position | Name           | Pre-Mixing Waiting Time (mm:ss) | Mixing Speed | Mixing Time (mm:ss) | Mixing Mode | Liquid Volume (μl) | Magnetic Attraction time (mm:ss) | Magnetic Attraction Count | Heating Temperature (°C) | Pause After Completion |
|------|---------------|----------------|---------------------------------|--------------|---------------------|-------------|--------------------|----------------------------------|---------------------------|--------------------------|------------------------|
| 1    | 6             | Transfer Beads | --                              | Fast         | 0:30                | Normal Mode | 930                | 0:30                             | 2                         | --                       | --                     |
| 2    | 2             | Store Beads    | --                              | Fast         | 0:30                | Normal Mode | 900                | 0:00                             | 0                         | --                       | --                     |
| 3    | 1             | Binding        | --                              | Fast         | 1:00                | Normal Mode | 900                | 0:00                             | 0                         | --                       | --                     |
| 4    | 2             | Transfer Beads | --                              | Fast         | 0:30                | Normal Mode | 930                | 0:30                             | 2                         | --                       | --                     |
| 5    | 1             | Binding        | --                              | Fast         | 3:00                | Normal Mode | 900                | 0:30                             | 0                         | --                       | --                     |
| 6    | 1             | Binding        | --                              | Fast         | 5:00                | Normal Mode | 900                | 0:30                             | 2                         | --                       | --                     |
| 7    | 2             | Wash 1         | --                              | Fast         | 0:30                | Normal Mode | 100                | 0:30                             | 0                         | --                       | --                     |
| 8    | 2             | Wash 2         | --                              | Fast         | 4:00                | Normal Mode | 900                | 0:30                             | 2                         | --                       | --                     |
| 9    | 3             | Wash 3         | --                              | Fast         | 0:30                | Normal Mode | 100                | 0:30                             | 0                         | --                       | --                     |
| 10   | 3             | Wash 4         | --                              | Fast         | 4:00                | Normal Mode | 900                | 0:30                             | 2                         | --                       | --                     |
| 11   | 4             | Wash 5         | --                              | Fast         | 3:00                | Normal Mode | 900                | 0:30                             | 2                         | --                       | --                     |
| 12   | 6             | Wash 6         | --                              | Fast         | 3:00                | Normal Mode | 900                | 0:30                             | 2                         | --                       | --                     |
| 13   | 5             | Elution        | 8:00                            | Fast         | 8:00                | Normal Mode | 120                | 0:30                             | 3                         | 75                       | --                     |
| 14   | 6             | Discard Beads  | --                              | Fast         | 0:30                | Normal Mode | 300                | 0                                | 0                         | --                       | --                     |

- After completion of the automated extraction procedure, collect DNA from well 5 of the single-sample cartridge and store under appropriate conditions.

## Appendix

### Related Products

#### 1. Instrument and Accessories

| Product name                      | Packing Size   | Cat.No      |
|-----------------------------------|----------------|-------------|
| TGuide S16 Nucleic Acid Extractor | 1 set          | OSE-S16-AM  |
| TGuide Smart Magnetic Tip Comb    | 200 pieces/box | OSE-TGA-S03 |
| TGuide Single Sample Tank Bracket | 5 pieces/box   | OSE-TGA-S32 |

#### 2. TGuide Smart Reagents Kits

| Product name                          | Cartridge | Cat.No    |
|---------------------------------------|-----------|-----------|
| TGuide Smart Viral DNA/RNA Kit        | 48 preps  | GDP604-DE |
| TGuide Smart Universal DNA Kit        | 48 preps  | GDP605-DE |
| TGuide Smart Envir-DNA Kit            | 48 preps  | GDP613-DE |
| TGuide Smart Hi-Plant DNA Kit         | 48 preps  | GDP617-DE |
| TGuide Smart DNA Purification Kit     | 48 preps  | GDP642-DE |
| TGuide Smart Micro Cell Total RNA Kit | 48 preps  | GDP670-DE |
| TGuide Smart Universal RNA Kit        | 48 preps  | GDP671-DE |
| TGuide Smart Hi-Plant RNA Kit         | 48 preps  | GDP672-DE |
| TGuide Smart Processed Food DNA Kit   | 48 preps  | GDP680-DE |
| TGuide Smart Viral DNA/RNA Kit        | 96 preps  | GDP604-E  |
| TGuide Smart Universal DNA Kit        | 96 preps  | GDP605-E  |
| TGuide Smart Envir-DNA Kit            | 96 preps  | GDP613-E  |
| TGuide Smart Hi-Plant DNA Kit         | 96 preps  | GDP617-E  |
| TGuide Smart DNA Purification Kit     | 96 preps  | GDP642-E  |
| TGuide Smart Micro Cell Total RNA Kit | 96 preps  | GDP670-E  |
| TGuide Smart Universal RNA Kit        | 96 preps  | GDP671-E  |
| TGuide Smart Hi-Plant RNA Kit         | 96 preps  | GDP672-E  |
| TGuide Smart Processed Food DNA Kit   | 96 preps  | GDP680-E  |